

AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT PCR kit



For Professional Use Only

Instruction Manual

KEY TO SYMBOLS USED

	Catalogue number		Caution
	Batch code		Contains sufficient for <n> tests
	In vitro diagnostic medical device		Use-by-date
	Version		Consult instructions for use
	Temperature limit		Keep away from sunlight
	Manufacturer	NCA	Negative control of amplification
	Date of manufacture	C-	Negative control of extraction
	Authorized representative in the European Community	C+Rotavirus / Astrovirus C+Norovirus genotype 2 / ST1	Positive controls of amplification
		IC	Internal control

1. INTENDED USE

AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT PCR kit is an *in vitro* nucleic acid amplification test for qualitative detection and differentiation of rotavirus A, norovirus genotype 2, and astrovirus RNA in the environmental samples (water sample concentrates) and clinical material (feces) using real-time hybridization-fluorescence detection of amplified products.

NOTE: The results of PCR analysis are taken into account in complex diagnostics of disease.

2. PRINCIPLE OF PCR DETECTION

Detection of rotavirus A, norovirus genotype 2, and astrovirus RNA includes RNA extraction from test samples and reverse transcription of RNA into cDNA combined with real-time PCR amplification of cDNA (RT-PCR). In the real-time PCR, the amplified product is detected with the use of fluorescent dyes. These dyes are linked to oligonucleotide probes, which bind specifically to the amplified product during thermocycling. The real-time monitoring of fluorescence intensities during the real-time PCR allows the detection of accumulating product without re-opening the reaction tubes after the PCR run.

AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT PCR kit is a qualitative test that contains the Internal Control (Internal Control STI-87-rec (IC)). It must be used in the extraction procedure in order to control the extraction process of each individual sample and to identify possible reaction inhibition.

AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT PCR kit uses "hot-start," which greatly reduces the frequency of nonspecifically primed reactions. "Hot-start" is guaranteed by separation of nucleotides and Taq-polymerase by a chemically modified polymerase (TaqF). The chemically modified polymerase (TaqF) is activated by heating at 95 °C for 15 min.

The results of amplification are registered in the following fluorescence channels:

Table 1

Channel for fluorophore	FAM	JOE
DNA-target		
RT-PCR-mix-1-FEP/FRT Rotavirus / Astrovirus	Rotavirus grA RNA	Astrovirus RNA
RT-PCR-mix-1-FEP/FRT Norovirus / ST1	Internal Control cDNA	Norovirus G2 RNA
Target gene		
RT-PCR-mix-1-FEP/FRT Rotavirus / Astrovirus	NSP2	Gene for capsid protein
RT-PCR-mix-1-FEP/FRT Norovirus / ST1	Artificially synthesized sequence	Gene for capsid protein

3. CONTENT

AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT PCR kit is produced in 1 form: variant FRT-50 F, R-V40(RG,iQ)-CE.

Variant FRT-50 F includes:

Reagent	Description	Volume, ml	Quantity
RT-PCR-mix-1-FEP/FRT Rotavirus / Astrovirus	clear liquid from colorless to light lilac colour	0.6	1 tube
RT-PCR-mix-1-FEP/FRT Norovirus / ST1	clear liquid from colorless to light lilac colour	0.6	1 tube
RT-PCR-mix-2-FEP/FRT	colorless clear liquid	0.3	2 tubes
Polymerase (TaqF)	colorless clear liquid	0.03	2 tubes
TM-Revertase (MMLv)	colorless clear liquid	0.015	2 tubes
RT-G-mix-2	colorless clear liquid	0.015	2 tubes
Positive Control cDNA Rotavirus-Flu / Astrovirus (C+ Rotavirus / Astrovirus)	colorless clear liquid	0.1	1 tube
Positive Control cDNA Norovirus genotype 2-Flu/ST1 (C+ Norovirus genotype 2 / ST1)	colorless clear liquid	0.1	1 tube
DNA-buffer	colorless clear liquid	0.5	1 tube
Negative Control (C-)*	colorless clear liquid	1.6	1 tube
Internal Control STI-87-rec (IC)**	colorless clear liquid	0.12	5 tubes
RNA-eluent***	colorless clear liquid	1.2	5 tubes

* must be used in the extraction procedure as Negative Control of Extraction

** add 10 µl of Internal Control STI-87-rec (IC) during the RNA extraction procedure directly to the sample/lysis mixture (see the RIBO-sorb K2-1-Et-50-CE, or RIBO-prep K2-9-Et-50-CE protocols).

***must be used in the extraction procedure.

Variant FRT-50 F is intended for 55 reactions (including controls).

4. ADDITIONAL REQUIREMENTS

- RNA extraction kit.
- Disposable powder-free gloves and a laboratory coat.
- Pipettes (adjustable).
- Sterile RNase-free pipette tips with aerosol filters (up to 100 µl).
- Tube racks.
- Vortex mixer.
- Desktop centrifuge with rotor for 2 ml reaction tubes.
- PCR box.
- Real-time instruments (for example, Rotor-Gene 3000/6000 (Corbett Research, Australia), Rotor-Gene Q (QIAGEN, Germany), iCycler iQ5 (Bio-Rad, USA), Mx3000P (Stratagene, USA)).
- Disposable polypropylene PCR tubes:
 - a) 0.2-ml tube (flat cap, nonstriped) for 36-well rotor if a rotor-type instrument is used.
 - b) 0.2-ml tube (domed cap) if plate-type instrument is used.
- Refrigerator for 2–8 °C.
- Deep-freezer at the temperature from minus 24 to minus 16 °C.
- Reservoir for used tips.

5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol filters and use a new tip for every procedure.
- Store all extracted positive material (specimens, controls and amplicons) away from all other reagents and add it to the reaction mix in a distantly separated facility.
- Thaw all components thoroughly at room temperature before starting an assay.
- When thawed, mix the components and centrifuge briefly.
- Use disposable protective gloves and laboratory cloths, and protect eyes while samples and reagents handling. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local regulations.
- Samples should be considered potentially infectious and handled in biological cabinet in compliance with appropriate biosafety practices.
- Clean and disinfect all samples or reagents spills using a disinfectant, such as 0.5 % sodium hypochlorite or another suitable disinfectant.
- Avoid inhalation of vapors, samples and reagents contact with the skin, eyes, and mucous membranes. Harmful if swallowed. If these solutions come into contact, rinse the injured area immediately with water and seek medical advice if necessary.
- Safety Data Sheets (SDS) are available on request.
- Use of this product should be limited to personnel trained in DNA amplification techniques.
- Workflow in the laboratory must be one-directional, beginning in the Extraction Area and moving to the Amplification and Detection Area. Do not return samples, equipment and reagents in the area where the previous step was performed.



Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

6. SAMPLING AND HANDLING

Obtaining samples of biological materials for PCR-analysis, transportation, and storage are described in the manufacturer's handbook [1]. It is recommended that this handbook is read before starting work.

AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT PCR kit is intended for analysis of the RNA extracted with RNA extraction kits from:

- water sample concentrates (pretreatment is not required),
- feces (pretreatment should be carried out as described in manufacturer's handbook [1]).

7. WORKING CONDITIONS

AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT PCR kit should be used at 18–25 °C.

8. PROTOCOL

8.1. RNA extraction

It is recommended that the following nucleic acid extraction kits are used:

- RIBO-prep **REF** K2-9-Et-50-CE;
- RIBO-sorb **REF** K2-1-Et-50-CE.

The RNA extraction of each test sample is carried out in the presence of **Internal Control STI-87-rec (IC)**.

In the extraction procedure it is necessary to carry out the control reaction as follows:

- **C-** – Add **100 µl of Negative Control (C-)** to the tube labelled **C-** (Negative Control of Extraction).

NOTE: Extract RNA according to the manufacturer's protocols.

NOTE: In case of extracting with RIBO-sorb reagent kit, the volume of **Internal Control STI-87-rec (IC)** reagent added to each tube is **10 µl**.

NOTE: Use RNA-eluent included in this PCR kit during RNA extraction.

8.2. Preparing RT-PCR

8.2.1. Preparing tubes for RT-PCR

The total reaction volume is **25 µl**, the volume of cDNA sample is **10 µl**.

The type of tubes depends on the type of PCR real-time instrument.

Use disposable filter tips for adding reagents, cDNA and control samples into tubes.

1. Reaction mixture components should be mixed just before analysis with calculating for the required reaction number (including test and control samples) according to Table 2. Note that even for analysis of one test or control RNA sample it is necessary to carry out all controls of RT-PCR (positive (C+) and negative (NCA)) for each RT-PCR-mix-1. It is recommended to mix the reagents for an even reaction number to ensure more exact dosage.
2. Take the required number of tubes for amplification of test and control samples.
3. To prepare the reaction mixture, mix one of the **RT-PCR-mix-1-FEP/FRT (RT-PCR-mix-1-FEP/FRT Rotavirus / Astrovirus or RT-PCR-mix-1-FEP/FRT Norovirus / STI)**, **RT-PCR-mix-2 FEP/FRT, polymerase (TaqF)**, and **TM-Revertase (MMIv)** according to Table 1. Vortex the tubes thoroughly. Make sure that there are no drops on the walls of the tubes.
4. Transfer **15 µl** of the prepared mixture to the prepared tubes. Dispose of the unused reaction mixture.

Table 2

Scheme of reaction mixture preparation

Reagent volume per one reaction, µl		Reagent volume for specified number of reactions, µl				
		10.00	5.00	0.25	0.50	0.25
Number of test samples	Number of reactions ¹	RT-PCR-mix-1-FEP/FRT	RT-PCR-mix-2-FEP/FRT	RT-G-mix-2	Polymerase (TaqF)	TM-Revertase (MMIv)
2	6	60	30	1.5	3.0	1.5
4	8	80	40	2.0	4.0	2.0
6	10	100	50	2.5	5.0	2.5
8	12	120	60	3.0	6.0	3.0
10	14	140	70	3.5	7.0	3.5
12	16	160	80	4.0	8.0	4.0
14	18	180	90	4.5	9.0	4.5
16	20	200	100	5.0	10.0	5.0
18	22	220	110	5.5	11.0	5.5
20	24	240	120	6.0	12.0	6.0
22	26	260	130	6.5	13.0	6.5
24	28	280	140	7.0	14.0	7.0
26	30	300	150	7.5	15.0	7.5
28	32	320	160	8.0	16.0	8.0

5. Add **10 µl** of RNA obtained at the RNA extraction stage to the prepared tubes using tips with aerosol filter.

6. Carry out the control amplification reactions:

- **NCA** – Add **10 µl of DNA-buffer** to the tube labeled **NCA** (Negative Control of Amplification).
- **C+Rotavirus / Astrovirus** – Add **10 µl of Positive Control cDNA Rotavirus-Flu / Astrovirus (C+Rotavirus / Astrovirus)** (in case of using RT-PCR-mix-1-FEP/FRT **Rotavirus / Astrovirus**) to the tube labeled **C+Rotavirus / Astrovirus** (Positive Control of Amplification)
- **C+Norovirus genotype 2 / STI** – Add **10 µl of Positive Control cDNA Norovirus genotype 2-Flu / STI (C+Norovirus genotype 2 / STI)** (in case of using RT-PCR-mix-1-FEP/FRT **Norovirus / STI**) to the tube labeled **C+Norovirus genotype 2 / STI** (Positive Control of Amplification)
- **C-** – Add **10 µl of the sample extracted from the Negative Control (C-)** reagent to the tube labeled **C-** (Negative control of Extraction).

NOTE: Avoid transferring sorbent beads together with the RNA sample in case of extraction with the RIBO-sorb kit.

¹ Test samples (N) + control of RNA extraction + 2 controls of RT-PCR + extra reaction (N+1+2+1).

8.2.2. Amplification

1. Create a temperature profile on your instrument as follows:

Table 3

Step	Rotor-type Instruments ²			Plate-type Instruments ³		
	Temperature, °C	Time	Cycles	Temperature, °C	Time	Cycles
1	50	30 min	1	50	30 min	1
2	95	15 min	1	95	15 min	1
3	95	10 s	45	95	10 s	45
	60	25 s Fluorescence acquiring		60	25 s Fluorescence acquiring	
	72	10 s		72	10 s	

Fluorescent signal is detected in the channels for the FAM and JOE fluorophores (other channels are enabled if several tests are simultaneously carried out in a single run).

2. Adjust the fluorescence channel sensitivity according to *Important Product Information Bulletin* and Guidelines [2].
3. Insert tubes into the reaction module of the device.
4. Run the amplification program with fluorescence detection.
5. Analyze results after the amplification program is completed.

9. DATA ANALYSIS

Analysis of results is performed by the software of the real-time PCR instrument used by measuring fluorescence signal accumulation in two channels (see table 4).

Table 4

Correspondence table of detection channels, RT-PCR-mixes-1 and pathogens

Channel for fluorophore	RT-PCR-mix-1-FEP/FRT Rotavirus / Astrovirus	RT-PCR-mix-1-FEP/FRT Norovirus / STI
FAM	Rotavirus grA cDNA	Internal Control STI-87-rec
JOE	Astrovirus cDNA	Norovirus G2 cDNA

Result interpretation

Results are interpreted by the crossing (or not-crossing) the fluorescence curve with the threshold line set at the specific level that corresponds to the presence (or absence) of a Ct value of the RNA sample in the corresponding column of the results grid.

Results should be interpreted in accordance with Table 5 and *Important Product Information Bulletin*.

Table 5

Interpretation of results

Channel for fluorophore	RT-PCR-mix-1-FEP/FRT Rotavirus / Astrovirus	RT-PCR-mix-1-FEP/FRT Norovirus / STI
FAM	< boundary value Rotavirus grA RNA is detected	< boundary value IC cDNA is detected. The result of sample is valid
	Absent or > boundary value Rotavirus grA RNA is not detected ⁴	Absent or > boundary value Invalid result ⁵
JOE	< boundary value Astrovirus RNA is detected	< boundary value Norovirus G2 RNA is detected
	Absent or > boundary value Astrovirus RNA is not detected ⁴	Absent or > boundary value Norovirus G2 RNA is not detected ⁴

NOTE: Boundary Ct values are specified in the *Important Product Information Bulletin* enclosed to the PCR kit.

The result of the analysis is considered reliable only if the results obtained for the Positive and Negative Controls of amplification as well as for the Negative Control of extraction are correct (see Table 6).

Table 6

Results for controls

RT-PCR-mix-1	Control	Stage for control	Ct value in the channel for fluorophore	
			FAM	JOE
RT-PCR-mix-1-FEP/FRT Norovirus / STI	C-	RNA extraction	≤ boundary value	Absent or > boundary value
	NCA	PCR	Absent or > boundary value	Absent or > boundary value
	C+Norovirus genotype 2 / STI	PCR	< boundary value	< boundary value
RT-PCR-mix-1-FEP/FRT Rotavirus / Astrovirus	C-	RNA extraction	Absent or > boundary value	Absent or > boundary value
	NCA	PCR	Absent or > boundary value	Absent or > boundary value
	C+Rotavirus / Astrovirus	PCR	< boundary value	< boundary value

10. TROUBLESHOOTING

Results of analysis are not taken into account in the following cases:

1. If the Ct value determined for the Positive Control of amplification (C+) in the channels for the JOE and FAM fluorophores is greater than the boundary Ct value, the amplification and detection should be repeated for all samples in which the signal in the channels for the JOE and FAM fluorophores was greater than the boundary value with the appropriate RT-PCR-mix.
2. If the signal for the Negative control of extraction (C-) (except for RT-PCR-mix-1-FEP/FRT **Norovirus / STI**) and/or the Negative control of amplification (NCA) in the channels for the JOE and FAM fluorophores is less than the boundary value, PCR should be repeated (starting from RNA extraction stage) for all samples in which the pathogen cDNA was detected.

If you have any further questions or if you encounter problems, please contact our Authorized representative in the European Community.

² Rotor-Gene 3000, Rotor-Gene 6000, or equivalent.

³ iCycler iQ5, Mx3000P, or equivalent.

⁴ Only if the Ct value for RT-PCR-mix-FEP/FRT **Norovirus / STI** in the FAM channel is less than the boundary value.

⁵ If Ct value for RT-PCR-mix-FEP/FRT **Norovirus / STI** in the FAM channel is absent or greater than the boundary value, the negative result obtained with the other PCR-mix-1 is considered invalid; therefore, the sample should be examined once again starting from RNA extraction.

11. TRANSPORTATION

AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT PCR kit should be transported at 2–8 °C for no longer than 5 days.

12. STABILITY AND STORAGE

All components of the **AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT** PCR kit are to be stored at 2–8 °C when not in use (except for RT-PCR-mix-1-FEP/FRT *Rotavirus / Astrovirus*, RT-PCR-mix-1-FEP/FRT *Norovirus / STI*, RT-PCR-mix-2-FEP/FRT, polymerase (TaqF), TM-Revertase (MMIv), and RT-G-mix-2) are to be stored at temperature from minus 24 to minus 16 °C. The shelf life of opened reagents is the same as that of unopened reagents, unless otherwise stated.

NOTE: RT-PCR-mix-1-FEP/FRT *Rotavirus / Astrovirus*, RT-PCR-mix-1-FEP/FRT *Norovirus / STI*, RT-PCR-mix-2-FEP/FRT, polymerase (TaqF), TM-Revertase (MMIv), and RT-G-mix-2 are to be stored at temperature from minus 24 to minus 16 °C.

NOTE: RT-PCR-mix-1-FEP/FRT *Rotavirus / Astrovirus* and RT-PCR-mix-1-FEP/FRT *Norovirus / STI* are to be kept away from light

13. SPECIFICATIONS

13.1. Sensitivity

The analytical sensitivity of **AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT** PCR kit is specified in the table below.

Pathogen	Test material	RNA/DNA extraction kit	PCR kit	Analytical sensitivity, GE/ml ⁶
<i>Rotavirus A</i>	Feces	RIBO-prep	PCR kit variant FRT-50 F	1 x 10 ⁴
<i>Norovirus genotype 2</i>	Feces	RIBO-prep	PCR kit variant FRT-50 F	5 x 10 ³
<i>Astrovirus</i>	Feces	RIBO-prep	PCR kit variant FRT-50 F	1 x 10 ⁴

13.2. Specificity

The analytical specificity of **AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT** PCR kit is ensured by selection of specific primers and probes as well as strict reaction conditions. The primers and probes have been checked for possible homologies to all sequences published in gene banks by sequence comparison analysis.

Specificity was confirmed on the following microorganism strains: *Enterovirus* strains (Coxsackie B1, B2, B3, B4, B5, and B6; Polio (Sabin) I, II, and III); *Adenovirus* serogroups 5 and 7; *influenza virus A* (H13N2, H9N2, H8N4, H2N3, H4N6, H11N6, H12N5, H3N8, H1N1, H6N2, H10N7, and H5N1) and B; *rhinoviruses*; *RS viruses*; human *adenovirus* types 3, 5, 7, 37, and 40; *Salmonella enteritidis* S-6, *S. choleraesuis* 370, *S. typhimurium* 371, *S. dublin* 373, *S. typhi* C1, *S. abortusovis* 372, and *S. gallinarum-pullorum*; *Shigella flexneri* 851b; *Campylobacter fetus* ssp. *fetus* 25936 and *C. jejuni* ssp. *jejuni* 43435; *Clebsiella K 65* SW4; *Listeria monocitogenes* USKHCH 19 and *L. monocitogenes* USKHCH 52; *Proteus vulgaris* 115/98; *Pseudomonas aeruginosa* DN c1; *Staphylococcus aureus* 653 and *S. aureus* 29112; *Morganella morgani* 619 c 01; and *Enterobacter faecalis* 356; as well as 44 *Norovirus* isolates of different gene clusters of genotypes 1 and 2; 40 *Rotavirus* strains of different [PJG] types, 19 *Astrovirus* strains of serogroups 1, 2, 4, 5, and 8; and 15 *Adenovirus* strains of different types and the following bacterial strains (see table 7).

Table 7

The panel of bacterial pathogens
Center for Disease Control and Prevention (CDC, USA)

Strain ID	Organism	Strain ID	Organism
K2033	<i>Salmonella ser. Grumpensis</i>	K2015	<i>Salmonella ser. Oranienburg</i>
K1806	<i>Salmonella ser. Newport</i>	AM01144	<i>Salmonella ser. Newport</i>
K2077	<i>Salmonella ser. Enteritidis</i>	K1810	<i>Salmonella ser. Anatum</i>
83-99	<i>Salmonella ser. Typhimurium</i>	K1991	<i>Salmonella ser. Typhimurium</i>
PS505	<i>Shigella boydii</i>	K1898	<i>Salmonella ser. Heidelberg</i>
PS408	<i>Shigella sonnei</i>	PS555	<i>Shigella boydii</i>
B4003	<i>Shigella sonnei</i>	F6446	<i>Shigella dysenteriae</i>
PS801	<i>Shigella dysenteriae</i>	S821X1	<i>Shigella dysenteriae type 1</i>
C898	<i>Shigella dysenteriae type 1</i>	S177X1	<i>Shigella dysenteriae type 1</i>
F2035	<i>Shigella flexneri</i>	S3314	<i>Shigella dysenteriae type 2</i>
E2539-C1	Enterotoxigenic <i>Escherichia coli</i> (ETEC)	PS071	<i>Shigella flexneri</i>
H10407	Enterotoxigenic <i>Escherichia coli</i> (ETEC)	PS050	<i>Shigella flexneri</i>
F1008	Enterotoxigenic <i>Escherichia coli</i> (ETEC)	F7862	<i>Shigella flexneri</i>
EDL 933	Shiga-toxin <i>E. coli</i> (STEC)	TX1	Enterotoxigenic <i>Escherichia coli</i> (ETEC)
3543-01	Shiga-toxin <i>E. coli</i> (STEC)	3525-01	Shiga-toxin <i>Escherichia coli</i> (STEC)
4752-71	<i>Proteus vulgaris</i>	25922	<i>Escherichia coli</i> O6:H1
QA/QC	<i>Citrobacter freundii</i>	621-64	<i>Citrobacter freundii</i>
QA/QC	<i>Aeromonas</i>	3910-68	<i>Aeromonas spp.</i>
3043-74	<i>Serratia marcescens</i>	E9113	<i>Vibrio cholerae</i>
QA/QC	<i>Serratia marcescens</i>	501-83	<i>Edwardsiella spp.</i>
F7894	<i>Vibrio vulnificus</i>	587-82	<i>Providencia stuartii</i>
F8515	<i>Yersinia enterocolitica</i>	27853	<i>Pseudomonas aeruginosa</i>
F8510	<i>Yersinia enterocolitica</i>	D4989	<i>Helicobacter cinaedi</i>
K4299	<i>Vibrio parahaemolyticus</i>	D6827	<i>Helicobacter pullorum</i>
F9835	<i>Vibrio parahaemolyticus</i>	D5127	<i>Helicobacter pylori</i>
K2023	<i>Salmonella ser. Kentucky</i>	D2686	<i>Arcobacter butzleri</i>
K1684	<i>Salmonella O-1, 4, 12 gr. B</i>		

There were no nonspecific test responses during examination of human DNA as well as a DNA panel of the above-mentioned microorganisms.

The clinical specificity of **AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT** PCR kit was confirmed in laboratory clinical trials.

⁶ Genome equivalents (GE) of the microorganism per 1 ml of a sample.

14. REFERENCES

- Handbook "Sampling, Transportation, and Storage of Clinical Material for PCR Diagnostics", developed by Federal State Institute of Science "Central Research Institute of Epidemiology" of Federal Service for Surveillance on Consumers' Rights Protection and Human Well-Being.
- Guidelines to the **AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT** PCR kit for qualitative detection and differentiation of *Rotavirus A*, *Norovirus genotype 2*, and *Astrovirus* RNA in environmental samples and clinical material by the polymerase chain reaction (PCR) with real-time hybridization-fluorescence detection developed by Federal Budget Institute of Science "Central Research Institute for Epidemiology".

15. QUALITY CONTROL

In accordance with Federal Budget Institution of Science "Central Research Institute for Epidemiology" ISO 13485-Certified Quality Management System, each lot of **AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT** PCR kit has been tested against predetermined specifications to ensure consistent product quality.

List of Changes Made in the Instruction Manual

VER	Location of changes	Essence of changes
25.06.11 LA	Cover page, text	The name of Institution was changed to Federal Budget Institution of Science "Central Research Institute for Epidemiology"
27.04.15 ME	Through the text	Corrections according to the template. Grammar corrections
	8.1. RNA Extraction	Information about controls of extraction was added
	8.2.1 Preparing tubes for RT-PCR	Appendix 1 was integrated into the text of the instruction manual as Table 1
	10. Troubleshooting	The section was rewritten
21.03.18 PM	14. References	The reference to Guidelines was added
	3. Content	The color of reagents was specified
19.05.20 EM	Through the text	The text formatting was changed
	Footer	The phrase "Not for use in the Russian Federation" was added
11.03.21 MM	2. Principle of PCR detection	The table with targets was added
	—	The name, address and contact information for Authorized representative in the European Community was changed

AmpliSens®



Ecoli Dx, s.r.o., Purkyňova 74/2
110 00 Praha 1, Czech Republic
Tel.: +420 325 209 912
Cell: +420 739 802 523



Federal Budget Institute of Science "Central Research Institute for Epidemiology"
3A Novogireevskaya Street
Moscow 111123 Russia