



For Professional Use Only

# AmpliSens<sup>®</sup> *U.parvum* / *U.urealyticum*-FRT

## PCR kit

### Instruction Manual

# AmpliSens<sup>®</sup>



Ecoli s.r.o., Studenohorska 12  
841 03 Bratislava 47  
Slovak Republic  
Tel.: +421 2 6478 9336  
Fax: +421 2 6478 9040



Federal Budget Institute of  
Science "Central Research  
Institute for Epidemiology"  
3A Novogireevskaya Street  
Moscow 111123 Russia

## TABLE OF CONTENTS

1. INTENDED USE.....	3
2. PRINCIPLE OF PCR DETECTION .....	3
3. CONTENT.....	3
4. ADDITIONAL REQUIREMENTS .....	4
5. GENERAL PRECAUTIONS .....	5
6. SAMPLING AND HANDLING .....	6
7. WORKING CONDITIONS .....	6
8. PROTOCOL.....	6
9. DATA ANALYSIS .....	8
10. TROUBLESHOOTING .....	9
11. TRANSPORTATION .....	10
12. STABILITY AND STORAGE .....	10
13. SPECIFICATIONS .....	10
14. REFERENCES.....	11
15. QUALITY CONTROL .....	11
16. KEY TO SYMBOLS USED .....	12

## 1. INTENDED USE

**AmpliSens<sup>®</sup> *U.parvum* / *U.urealyticum*-FRT** PCR kit is an *in vitro* nucleic acid amplification test for qualitative detection and differentiation of *Ureaplasma parvum* and *Ureaplasma urealyticum* DNA in the clinical material (urogenital, rectal, and oropharyngeal swabs; conjunctival discharge; urine samples, prostate gland secretion) using real-time hybridization-fluorescence detection of amplified products.



The results of PCR analysis are taken into account in complex diagnostics of disease.

## 2. PRINCIPLE OF PCR DETECTION

*U.parvum* / *U.urealyticum* detection by the polymerase chain reaction (PCR) is based on the amplification of the pathogen genome specific region using specific *U.parvum* / *U.urealyticum* primers. In the real-time PCR, the amplified product is detected with the use of fluorescent dyes. These dyes are linked to oligonucleotide probes which bind specifically to the amplified product during thermocycling. The real-time monitoring of fluorescence intensities during the real-time PCR allows the detection of accumulating product without re-opening the reaction tubes after the PCR run.

**AmpliSens<sup>®</sup> *U.parvum* / *U.urealyticum*-FRT** PCR kit is a qualitative test that contains the Internal Control (Internal Control-FL (IC)). It must be used in the extraction procedure in order to control the extraction process of each individual sample and to identify possible reaction inhibition.

**AmpliSens<sup>®</sup> *U.parvum* / *U.urealyticum*-FRT** PCR kit uses “hot-start,” which greatly reduces the frequency of nonspecifically primed reactions. In variant FRT, “hot-start” is guaranteed by the separation of nucleotides and Taq-polymerase using a wax layer. Wax melts and reaction components mix only at 95 °C. In variant FRT-100 F, “hot-start” is guaranteed by the separation of nucleotides and Taq-polymerase using chemically modified polymerase (TaqF). The chemically modified polymerase (TaqF) is activated by heating at 95 °C for 15 min.

## 3. CONTENT

**AmpliSens<sup>®</sup> *U.parvum* / *U.urealyticum*-FRT** PCR kit is produced in 2 forms:

**AmpliSens<sup>®</sup> *U.parvum* / *U.urealyticum*-FRT** PCR kit variant FRT **REF** R-B19(RG)-CE.

**AmpliSens<sup>®</sup> *U.parvum* / *U.urealyticum*-FRT** PCR kit variant FRT-100 F, **REF** R-B19-F(RG,iQ)-CE.

**AmpliSens<sup>®</sup> *U.parvum* / *U.urealyticum*-FRT** PCR kit variant FRT includes:

<b>Reagent</b>	<b>Description</b>	<b>Volume, ml</b>	<b>Quantity</b>
<b>PCR-mix-1-FL <i>U.parvum</i> / <i>U.urealyticum</i></b> (ready-to-use single-dose test tubes ( <i>under wax</i> ))	clear liquid from colorless to light lilac colour	0.01	110 tubes of 0.2 ml
<b>PCR-mix-2-FL-red</b>	red clear liquid	1.1	1 tube
<b>Positive Control complex (C+)</b>	colorless clear liquid	0.2	1 tube
<b>DNA-buffer</b>	colorless clear liquid	0.5	1 tube
<b>Negative Control (C-)*</b>	colorless clear liquid	1.2	1 tube
<b>Internal Control-FL (IC)**</b>	colorless clear liquid	1.0	1 tube

\* must be used in the extraction procedure as Negative Control of Extraction.

\*\* add 10 µl of **Internal Control-FL (IC)** during the DNA extraction procedure directly to the sample/lysis mixture (see **DNA-sorb-AM** **REF** K1-12-100-CE protocol).

**AmpliSens® *U.parvum* / *U.urealyticum*-FRT** PCR kit is intended for 110 reactions (including controls).

**AmpliSens® *U.parvum* / *U.urealyticum*-FRT** PCR kit variant FRT-100 F includes:

<b>Reagent</b>	<b>Description</b>	<b>Volume, ml</b>	<b>Quantity</b>
<b>PCR-mix-1-FL <i>U.parvum</i> / <i>U.urealyticum</i></b>	clear liquid from colorless to light lilac colour	1.2	1 tube
<b>PCR-mix-2-FRT</b>	colorless clear liquid	0.3	2 tubes
<b>Polymerase (TaqF)</b>	colorless clear liquid	0.03	2 tubes
<b>Positive Control complex (C+)</b>	colorless clear liquid	0.2	1 tube
<b>DNA-buffer</b>	colorless clear liquid	0.5	1 tube
<b>Negative Control (C-)*</b>	colorless clear liquid	1.2	1 tube
<b>Internal Control-FL (IC)**</b>	colorless clear liquid	1.0	1 tube

\* must be used in the extraction procedure as Negative Control of extraction.

\*\* add 10 µl of **Internal Control-FL (IC)** during the DNA extraction procedure directly to the sample/lysis mixture (see **DNA-sorb-AM** **REF** K1-12-100-CE protocol).

**AmpliSens® *U.parvum* / *U.urealyticum*-FRT** PCR kit is intended for 110 reactions (including controls).

#### 4. ADDITIONAL REQUIREMENTS

- DNA extraction kit.

- Transport medium.
- Disposable powder-free gloves and a laboratory coat.
- Pipettes (adjustable).
- Sterile pipette tips with aerosol barriers (up to 100 µl).
- Tube racks.
- Vortex mixer.
- Desktop centrifuge with rotor for 2-ml reaction tubes.
- PCR box.
- Real-time instruments (for example, Rotor-Gene 3000/6000 (Corbett Research, Australia); Rotor-Gene Q (QIAGEN, Germany); iCycler iQ5 (Bio-Rad, USA); Mx3000P (Stratagene, USA).
- Disposable polypropylene PCR tubes (0.1- or 0.2-ml) when working with PCR kit variant FRT-100 F:
  - a) 0.2-ml PCR tubes with optical transparent domed or flat caps if a plate-type instrument is used;
  - b) 0.2-ml PCR tubes with flat caps or strips of four 0.1-ml Rotor-Gene PCR tubes if a rotor-type instrument is used.
- Refrigerator for 2–8 °C.
- Deep-freezer at the temperature from minus 24 to minus 16 °C.
- Reservoir for used tips.

## 5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol filters and use a new tip for every procedure.
- Store all extracted positive material (specimens, controls and amplicons) away from all other reagents and add it to the reaction mix in a distantly separated facility.
- Thaw all components thoroughly at room temperature before starting an assay.
- When thawed, mix the components and centrifuge briefly.
- Use disposable protective gloves and laboratory cloths, and protect eyes while samples and reagents handling. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local regulations.
- Samples should be considered potentially infectious and handled in biological cabinet in

compliance with appropriate biosafety practices.

- Clean and disinfect all samples or reagents spills using a disinfectant, such as 0.5 % sodium hypochlorite or another suitable disinfectant.
- Avoid samples and reagents contact with the skin, eyes, and mucous membranes. If these solutions come into contact, rinse the injured area immediately with water and seek medical advice immediately.
- Safety Data Sheets (SDS) are available on request.
- Use of this product should be limited to personnel trained in DNA amplification techniques.
- Workflow in the laboratory must be one-directional, beginning in the Extraction Area and moving to the Amplification and Detection Area. Do not return samples, equipment and reagents in the area where the previous step was performed.



Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

## 6. SAMPLING AND HANDLING



Obtaining samples of biological materials for PCR-analysis, transportation and storage are described in the manufacturer's handbook [1]. It is recommended that this handbook is read before starting work.

**AmpliSens<sup>®</sup> *U.parvum* / *U.urealyticum*-FRT** PCR kit is intended for analysis of the DNA extracted with DNA extraction kits from the clinical material (urogenital, rectal, and oropharyngeal swabs, conjunctival discharge, urine samples (a sediment of the first portion of the morning specimen), prostate gland secretion).

## 7. WORKING CONDITIONS

**AmpliSens<sup>®</sup> *U.parvum* / *U.urealyticum*-FRT** PCR kit should be used at 18–25 °C.

## 8. PROTOCOL

### 8.1. DNA extraction

It is recommended to use the following nucleic acid extraction kits:

- DNA-sorb-AM, **REF** K1-12-100-CE.
- For other nucleic acid extraction kits see Guidelines [2].

The DNA extraction of each test sample is carried out in the presence of **Internal Control-FL (IC)**.

In the extraction procedure it is necessary to carry out the control reactions as follows:

**C-** – Add **100 µl of Negative Control (C-)** to the tube labeled C-.



Extract DNA according to the manufacturer's protocol.

## 8.2. Preparing PCR

### 8.2.1. Preparing tubes for PCR

The type of tubes depends on the PCR instrument used for analysis. Use disposable filter tips for adding reagents, DNA and control samples into tubes.

#### Variant FRT

The total reaction volume is **30 µl**, the volume of DNA sample is **10 µl**.

1. Prepare the required number of tubes with **PCR-mix-1-FL *U.parvum* / *U.urealyticum*** and wax for amplification of DNA from clinical and control samples.
2. Add **10 µl of PCR-mix-2-FL-red** to the surface of the wax layer into each tube ensuring that it does not fall under the wax and mix with **PCR-mix-1-FL *U.parvum* / *U.urealyticum***.

#### Variant FRT-100 F

The total reaction volume is **25 µl**, the volume of DNA sample is **10 µl**.

1. Thaw the tube with **PCR-mix-2-FRT**. Vortex the tubes with **PCR-mix-1-FL *U.parvum* / *U.urealyticum*, PCR-mix-2-FRT, and polymerase (TaqF)** then centrifuge briefly.

Take the required number of the tubes/strips for amplification of DNA obtained from clinical and control samples.

2. For N reactions (including 2 controls of amplification) add to a new tube:

**10\*(N+1) µl of PCR-mix-1-FL *U.parvum* / *U.urealyticum*,**

**5.0\*(N+1) µl of PCR-mix-2-FRT,**

**0.5\*(N+1) µl of polymerase (TaqF).**

Vortex the tube, then centrifuge briefly. Transfer **15 µl** of the prepared mixture into each tube.

Steps 3 and 4 are carried out in both variants.

3. Add **10 µl of DNA samples** obtained at the DNA extraction stage.

4. Carry out the control amplification reactions:

**NCA** – Add **10 µl of DNA-buffer** to the tube labeled NCA (Negative Control of Amplification).

**C+** – Add **10 µl of Positive Control complex (C+)** to the tube labeled C+ (Positive control of amplification).

**C-** – Add **10 µl of the sample, extracted from Negative Control (C-) reagent** to the tube labeled C- (Negative Control of Extraction).

## 8.2.2. Amplification

1. Create a temperature profile on your instrument as follows:

Table 1

**AmpliSens-1 amplification program**

Step	Rotor-type Instruments <sup>1</sup>			Plate-type Instruments <sup>2</sup>		
	Temperature, °C	Time	Cycles	Temperature, °C	Time	Cycles
1	95	15 min	1	95	15 min	1
2	95	5 s	5	95	5 s	5
	60	20 s		60	20 s	
	72	15 s		72	15 s	
3	95	5 s	40	95	5 s	40
	60	20 s <i>Fluorescence acquiring</i>		60	30 s <i>Fluorescence acquiring</i>	
	72	15 s		72	15 s	

Fluorescent signal is detected in the channels for the FAM, JOE, and ROX fluorophores (other channels are enabled if several tests are simultaneously carried out in a single run).

2. Adjust the fluorescence channel sensitivity according to the *Important Product Information Bulletin* and Guidelines [2].
3. Insert tubes into the reaction module of the device.
4. Run the amplification program with fluorescence detection.
5. Analyze results after the amplification program is completed.

## 9. DATA ANALYSIS

Analysis of results is performed by the software of the real-time PCR instrument used by measuring fluorescence signal accumulation in three channels:

- The signal of the *U.parvum* DNA amplification product is detected in the channel for the FAM fluorophore.
- The signal of the *U.urealyticum* DNA amplification product is detected in the channel for the JOE fluorophore.
- The signal of the IC DNA amplification product is detected in the channel for the ROX fluorophore.

Results are interpreted by the crossing (or not crossing) the fluorescence curve with the threshold line set at the specific level that corresponds to the presence (or absence) of Ct value of the DNA sample in the corresponding column of the result grid.

Principle of interpretation is the following:

<sup>1</sup> For example, Rotor-Gene 3000, Rotor-Gene 6000, Rotor-Gene Q or equivalent.

<sup>2</sup> For example, iCycler iQ5, Mx3000P, Mx3000, DT-96 or equivalent.



- *Ureaplasma parvum* DNA is **detected** in a sample if the *Ct* value is determined in the result grid in the channel for the FAM fluorophore. Moreover, the fluorescence curve of the sample should cross the threshold line in the area of typical exponential growth of fluorescence.
- *Ureaplasma urealyticum* DNA is **detected** in a sample if the *Ct* value is determined in the result grid in the channel for the JOE fluorophore. Moreover, the fluorescence curve of the sample should cross the threshold line in the area of typical exponential growth of fluorescence.
- *Ureaplasma parvum* and *Ureaplasma urealyticum* DNA are **not detected** in a sample if the *Ct* value is not determined (absent) in the result grid in the channels for the FAM and JOE fluorophores (the fluorescence curve does not cross the threshold line), whereas the *Ct* value in the channel for ROX fluorophore is less than the specified boundary *Ct* value.
- The result is invalid if *Ct* value is not determined (absent) in the channels for FAM, JOE and ROX fluorophores. In such cases, the PCR analysis should be repeated for such samples.



Boundary *Ct* values are specified in the *Important Product Information Bulletin* enclosed to the PCR kit. See also Guidelines [2]

**The result of the analysis is considered reliable only if the results obtained for the Positive and Negative Controls of amplification as well as for the Negative Control of extraction are correct (see Table 2).**

Table 2

#### Results for controls

Control	Stage for control	<i>Ct</i> value in channel for fluorophore	
		FAM, JOE	ROX
<b>C-</b>	DNA extraction	Absent	<boundary value
<b>NCA</b>	PCR	Absent	Absent
<b>C+</b>	PCR	<boundary value	<boundary value

## 10. TROUBLESHOOTING

Results of analysis are not taken into account in the following cases:

1. If the *Ct* value determined for the Positive Control of Amplification (C+) in the channels for the FAM and/or JOE fluorophores is greater than the boundary *Ct* value or absent, the amplification should be repeated for all samples in which the boundary *Ct* is absent in the channels for the FAM and/or JOE fluorophores .

2. If the *Ct* value is determined for the Negative Control of Amplification (NCA) and/or Negative Control of Extraction (C-) in the channels for the FAM and/or JOE fluorophores, the PCR analysis should be repeated for all samples in which the *Ct* value is determined in the channels for the FAM and/or JOE fluorophores.

If you have any further questions or if you encounter problems, please contact our Authorized Representative in the European Community.

## 11. TRANSPORTATION

**AmpliSens® *U.parvum* / *U.urealyticum*-FRT** PCR kit should be transported at 2–8 °C for no longer than 5 days.

## 12. STABILITY AND STORAGE

All components of the **AmpliSens® *U.parvum* / *U.urealyticum*-FRT** PCR kit (except for polymerase (TaqF) and PCR-mix-2-FRT) are to be stored at 2–8 °C when not in use. All components of the **AmpliSens® *U.parvum* / *U.urealyticum*-FRT** PCR kit are stable until labeled expiry date. The shelf life of opened reagents is the same as that of unopened reagents, unless otherwise stated.



Polymerase (TaqF) and PCR-mix-2-FRT are to be stored at the temperature from minus 24 to minus 16 °C when not in use.



PCR-mix-1-FL *U.parvum/U.urealyticum* is to be kept away from light.

## 13. SPECIFICATIONS

### 13.1. Sensitivity

Clinical material	Nucleic acid extraction kit	Microorganism	Analytical sensitivity, GE/ml <sup>3</sup>
Urogenital swabs <sup>4</sup>	DNA-sorb-AM	<i>Ureaplasma parvum</i>	10 <sup>3</sup>
		<i>Ureaplasma urealyticum</i>	10 <sup>3</sup>
Urine (pretreatment is required)	DNA-sorb-AM	<i>Ureaplasma parvum</i>	5x10 <sup>3</sup>
		<i>Ureaplasma urealyticum</i>	5x10 <sup>3</sup>

### 13.2. Specificity

The analytical specificity of **AmpliSens® *U.parvum* / *U.urealyticum*-FRT** PCR kit is ensured by selection of specific primers and probes as well as stringent reaction

<sup>3</sup> Genome equivalents (GE) of the pathogen agent per 1 ml of a sample placed in the transport medium.

<sup>4</sup> Urogenital swabs are to be placed into Transport medium for swabs (**REF** 956-CE, **REF** 987-CE) or Transport Medium with Mucolytic Agent (**REF** 952-CE, **REF** 953-CE).

conditions. The primers and probes were checked for possible homologies to all sequences deposited in gene banks by sequence comparison analysis.

Nonspecific responses were absent while testing human DNA samples as well as a DNA panel of the following microorganisms: *Gardnerella vaginalis*, *Lactobacillus spp.*, *Escherichia coli*, *Staphylococcus spp.*, *Streptococcus spp.*, *Candida albicans*, *Mycoplasma hominis*, *Ureaplasma urealyticum*, *Ureaplasma parvum*, *Mycoplasma genitalium*, *Chlamydia trachomatis*, *Neisseria spp.*, *Trichomonas vaginalis*, *Neisseria gonorrhoeae*, *Treponema pallidum*, *Toxoplasma gondii*, HSV 1 and 2, CMV, and HPV.

The clinical specificity of **AmpliSens® U.parvum / U.urealyticum-FRT** PCR kit was confirmed in laboratory clinical trials.














#### 14. REFERENCES

1. Handbook “Sampling, Transportation, and Storage of Clinical Material for PCR Diagnostics”, developed by Federal Budget Institute of Science “Central Research Institute for Epidemiology” of Federal Service for Surveillance on Consumers’ Rights Protection and Human Well-Being, Moscow, 2010.
2. Guidelines “Real-Time PCR Detection of STIs and Other Reproductive Tract Infections”, developed by Federal Budget Institute of Science “Central Research Institute for Epidemiology” of Federal Service for Surveillance on Consumers’ Rights Protection and Human Well-Being, Moscow”

#### 15. QUALITY CONTROL

In compliance with Federal Budget Institute of Science “Central Research Institute for Epidemiology” ISO 13485-Certified Quality Management System, each lot of **AmpliSens® U.parvum / U.urealyticum-FRT** PCR kit has been tested against predetermined specifications to ensure consistent product quality.

## 16. KEY TO SYMBOLS USED

	Catalogue number		Caution
	Batch code		Sufficient for
	<i>In vitro</i> diagnostic medical device		Expiration Date
	Version		Consult instructions for use
	Temperature limitation		Keep away from sunlight
	Manufacturer	<b>NCA</b>	Negative control of amplification
	Date of manufacture	<b>C-</b>	Negative control of extraction
	Authorised representative in the European Community	<b>C+</b>	Positive control of amplification
		<b>IC</b>	Internal control

### List of Changes Made in the Instruction Manual

VER	Location of changes	Essence of changes
27.06.11 RT	Cover page, text	The name of Institute was changed to Federal Budget Institute of Science "Central Research Institute for Epidemiology"
30.11.15 PM	Through the text	Corrections in accordance with the template
	8.1. DNA extraction	Information about controls of extraction was added
	9. Data analysis 10. Troubleshooting	The sections were rewritten
18.12.17 PM	3. Content	The colour of the reagent was specified
15.03.18 PM	Footer, 3. Content	<b>REF</b> R-B19(iQ)-CE was deleted