



# ECOLI

MADE IN  
THE CZECH  
REPUBLIC  
(EU)

## COROTEST VTM

VIROLOGICAL TRANSPORT MEDIUM

Sampling set designed for taking clinical samples, storage and transport to the laboratory.

### USAGE

- Contains polypropylene tube with transport medium
- Includes throat and nasopharyngeal swabs
- Suitable for section sampling
- **Suitable for SARS-CoV-2 virus, causing COVID-19 disease**



### SPECIFICATIONS

#### Ingredients:

- Bovine serum fetale BSF
- Serum albumin
- Antibiotics for suppressing growth of microflora
- PBS buffer (pH = 7.3 at 25°C)
- Storage in darkness at room temperature below 25°C



Ver.01\_08/2020

### REGISTERED BRANCH

Ecoli s.r.o., odštěpný závod, Purkyňova 74/2, 110 00 PRAHA 1, Czech Republic

Phone: +420 325 209 912, Cell: +420 739 802 523, +420 731 233 938, e-mail: [ecoli@ecoli.sk](mailto:ecoli@ecoli.sk), ID: 08761183, VAT ID: CZ685028328

### HEADQUARTERS

Ecoli s.r.o., Studenohorská 12, 841 03 BRATISLAVA, Slovakia

Phone: +420 325 209 912, Cell: +420 739 802 523, +420 731 233 938, e-mail: [ecoli@ecoli.sk](mailto:ecoli@ecoli.sk), ID: 35 815 906

