

MEDICA 2023 INVITATION

November 13 – 16, Düsseldorf/Germany

Introducing eSens, the truly European AmpliSens®



- New range of PCR kits
- **Original AmpliSens® chemistry**
(the same composition of the solutions, the same content of the kits)
- **Now, made in EU!**
- CE IVD certification



eQuantia 48

- Real-Time PCR cycler for 48 x 0,2 ml tubes
- 5-channel optics (FAM, HEX, ROX, Cy5 and Cy5.5)
with separate LED for each channel



Fully automated nucleic acid extraction platform

- Liquid handling instrument based on magnetic particle extraction technology
- Purifies 1 to 12 samples simultaneously per run
- Range of purification kits that contains all prefilled reagents and consumables
incl. Covid-19 Viral RNA/DNA extraction kit



**DISTRIBUTORS
WANTED**



Take your time and come to us to:

Hall 1, Booth no. D01

Contact us and arrange an appointment:
ecolidx@ecolidx.com

Yours Sincerely,
František Šmíd
International Sales Manager

MEDICA

